size:450*143mm 20240423

PBSBT_manual_low-res

FRONT

PBSBT

Smart Body Scale

bluetooth connectivity

Free Android & iOS app
LCD Display with Blue Backlight
8 Users Memory

Bluetooth

USER MANUAL

USER GUIDE

EN

BLUETOOTH BODY FAT SCALE

For first use

1. Battery Installation
Open the battery door, put into the battery with correct polarity.
2. App Installation
Scan the QR code on the first page or download „OKOK International“ app from App Store or Android Store.
3. Registration
Turn on the bluetooth, turn on the OKOK, follow the App steps to register the user profile.

Twelve health index checking

1. Default eight
body fat
body hydration
body muscle
body bone
methabolism
weight
BMI
visceral fat
2. Additional four index, release them or not from setting
protein
obesity degree
body age
weight without fat

Multi Users

Total 8 users share one user account in one device.
Auto indentify users according to different weight memory.

Target weight

Users can set their own target weight to make the health plan.

Attention

1. Please make sure good contact between bare feet and metal electrodes, otherwise it will display your weight only.
2. Ensure the scale is on a flat, even and firm surface. Avoid carpet or soft surface.
3. The surface of the scale will be slippery if wet. Please make sure both the glass and your feet are dry before use.
4. This scale is not water proof. Do not expose to moisture.
5. Clean the scale with a slightly damp cloth. DO not use any chemical/abrasive cleaning agents to clean the scale.
6. Please remove the batteries, if the scale is not to being used for a long time.
7. The appliance is intended solely for domestic use within the home.
8. Children should be supervised to ensure that they do not play with the appliance.

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PLATINET S.A. 30-798 Krakow, Stiiwaka 48, PL www.platinet.eu